



ADAPTATION ASSEMBLY INSTRUCTIONS



AmTryke
The therapeutic tricycle for people with disabilities

For questions or to place an order, please contact:

Amtryke, LLC.
(800) 838-1845 x114
amtryke@ambucs.org
www.ambucs.org

Vertical Handgrips

Part Number #30-23-0343

Fits AM-10, AM-12S, AM-12, and AM-16



Tools Needed: 14mm wrench

Note! If the Tryke is unassembled skip to Step 2

Installation Instructions

1) Remove the Hand Grips from the Hand Cranks

Using a 14mm wrench loosen and unscrew each hand grip, making sure to note which is right hand threaded and which is left hand threaded. They should be stamped "R" and "L" on the end of the thread.

2) Install Vertical Handgrips

Located on the flat surface of threaded brass portion is stamped either "L" for left or "R" for right. The right side of the bike has the chain guard on. Find the vertical hand grip that has the "R" stamped on the flat portion. Thread this into the right crank arm "same side chain guard is on" turning it to the right "clockwise" until snug. Use a 14mm wrench to finish tightening. Thread the left vertical grip into the left crank arm turning it counter clockwise. Use a 14mm wrench to finish tightening.

Revised: 6/20/2016

n:\amtryke\manuals and instructions\accessories\vehicle handgrips\vertical hand grips.doc