



# ADAPTATION ASSEMBLY INSTRUCTIONS



The therapeutic tricycle for people with disabilities

For questions or to place an order, please contact:

Amtryke, LLC.  
(800) 838-1845 x114  
amtryke@ambucs.org  
www.ambucs.org

## AmTryke® Wrist Wrap

XS = 7 inches and smaller - # 30-21-0160XS

S = 8 inches - # 30-21-0160S

MED = 10 inches - # 30-21-0160M

LG = 11 inches - # 30-21-0160LG




XL = 12 inches - # 30-21-0160XL

### Wrist Wrap Directions

1. Place hand on tryke hand grip/handle bar with thumb placement below grip and fingers above hand grip.
2. Place small T or L end of wrist wrap around the top side of the wrist and secure hook and loop.
3. Pull material snug over back of hand toward underside of wrist.
4. Secure hook and loop fastener of strap around wrist / back of hand.
5. Discomfort, constriction of blood flow or red marks should be avoided.

Measure from point **A** to **B** for correct size.

Xtra Small	Small	Medium	Large	Xtra Large
7"	8"	10"	11"	12"



Revised: 12/28/2015

n:\amtryke\manuals and instructions\accessories\wrist wraps\wrist