



## AmTryke JT-2300-USS Part #: 50-FC-2300-USS

Installation Instructions for JT2300-USS (Under Seat Steering)

NOTE: These instructions are intended as a quick reference guide. Please be sure to use the JT-2000 Assembly Guide as the main reference. It contains all the frame and rear deck instructions, since these two models share the same basic frame.

*Use the JT-2000 Assembly Instructions for installation* of:

- Wheels
- Rear deck
- Seat
- Small drive chain
- Front brake, caliper and adjustment
- Pedals
- Reflectors

1) Check for all parts and make sure you have received the JT-2300-USS correct parts.



Parts Bags



Hardware Bag



## Reflectors



- 2) Remove main frame from the box.



- 3) Install rear deck (shown in proper position).



- 4) You may prefer to install rear tires first – before attaching to main frame.
- 5) See JT-2000 instructions for rear deck and chain assembly instructions.




- 6) Your bike should look like this before handles are installed.

Revised: 11/3/2015

n:\amtryke\manuals and instructions\models\jt-2300 uss\jt-2300 uss assembly instructions.docx



- 7) Insert the handles into each side of the USS steering tubes. You can adjust and tighten later to your rider's comfort.



- 8) The skewer for the front wheel looks like this before installation. One (1) spring goes on each side. Be sure to orient them as shown.



- 9) Here is your bike with the seat frame installed.
- 10) See the JT-2000 instructions for more information

Questions or Concerns Regarding your Build? No problem! Please Contact Derek Shaw,  
AmTryke Technical and Customer Support Coordinator @ 1-800-838-1845 x114 or  
[dereks@ambucs.org](mailto:dereks@ambucs.org)