



Amtryke Therapeutic Tricycle Request Form

(Must be filled out completely by adult rider or parent/guardian)

Recipient's Name: _____ Age: _____ Date of Request: _____
Mailing Address: _____ Phone #: _____
City/State/Zip: _____ Email: _____
Diagnosis: _____

If Recipient is Under Age 18

Parent/Guardian Name: _____
If different from above
Mailing Address: _____ Phone #: _____
City/State/Zip: _____ Email: _____

Secondary Contact Name: _____ Phone #: _____

Treating Therapist's Name: _____
Phone #: _____ Email: _____

How did you hear about the Amtryke Therapeutic Tricycle? (Check all that apply)

Therapist Website AMBUCS Member Other: _____

Will you need financial assistance to obtain the tricycle? Yes No

If yes, how much can you pay? _____

Note: Amtryke therapeutic tricycles are distributed based on available funds and need. Individual placements of Amtryke therapeutic tricycles are at the discretion of the local chapter or parent organization.

Tell us about the recipient*:

*This information will be made public to help obtain funding. Please don't include information you don't want shared.

Including a photo of the recipient will help us obtain a sponsor to help you pay for the Amtryke more quickly. Digital images preferred but we also accept professionally printed glossy photos. No photocopies or folded images. By including a photo, you are giving consent for AMBUCS to use the image online and in print to help obtain a funding.

Signature: _____ Date: _____

Individuals will not be considered for placement until all three forms are returned to AMBUCS: this form filled out by the recipient/guardian, as well as the Assessment Form and Tryke Selection Form filled out by the therapist.

Amtryke Therapeutic Tricycle Waiver Form

(Must be filled out completely by adult rider or parent/guardian)

AMBUCS members nationwide are dedicated to creating opportunities for mobility and independence for people with disabilities by providing Amtryke therapeutic tricycles, offering educational scholarships to therapy students and performing various forms of community service.

Purpose: The Amtryke therapeutic tricycle was designed for people with disabilities. It creates a feeling of freedom, builds self-esteem, strengthens muscles and improves motor coordination and range of motion—all while making exercise fun.

Steering: Initially, the rider may have difficulty turning or changing directions. Encourage the rider to go straight ahead, back up and slowly turn around. On many models there are three steering options for the Amtryke. On the front column of the tricycle you will find two holes for the steering pin. The top hole is straight steering, the bottom hole allows a 20 degree turning radius. Leaving the pin out gives the rider free steering.

Safety Cautions

- Fast speeds and sharp turns can cause the Amtryke therapeutic tricycle to tip or turn over.
- Always wear a helmet when riding an Amtryke. Use of other protective gear is highly recommended.
- Adult supervision required if used by younger or developmentally delayed riders.
- Use caution near vehicles, swimming pools and other bodies of water, hills, alleys and sloped driveways.
- Always wear shoes.
- Never allow more than one rider.
- Use of the steering pin is recommended to prevent over-steering or possible tip-overs.

The information contained in this service is not intended nor implied by National AMBUCS™, Inc. to be professional medical advice by National AMBUCS, Inc. Always seek the advice of your physician, therapist or other qualified healthcare provider prior to starting any treatment or with any question you may have regarding a medical condition. Nothing contained in this document is intended by National AMBUCS, Inc. to be for medical diagnosis or treatment by National AMBUCS, Inc. or on behalf of National AMBUCS, Inc.

In no event shall National AMBUCS, Inc. be liable for any direct, indirect, incidental, consequential, special, exemplary, punitive, or any other monetary or other damages, personal injury or property damages, fees, fines, costs, attorney fees, or liabilities of any kind arising out of or relating in any way to this service or use of the Amtryke® therapeutic tricycle, and/or content or information provided herein.

I agree that the rider listed on the front of this document may be photographed. I also agree that the photos and recipient's name may be used in promotional efforts for National AMBUCS, Amtryke or the local AMBUCS chapter. I further grant AMBUCS the ability to use the photos and name for advertising/publicity purposes without additional compensation, except where prohibited by law. If anybody in my party does not want to be photographed under these same terms, I will let the photographer know as soon as possible.

By signing below, I acknowledge that I have read and understood this liability waiver.

Recipient's Name: _____

Adult Recipient Signature: _____

If Recipient is Under Age 18

Legal Guardian Name: _____

Legal Guardian Signature: _____ Date: _____

Individuals will not be considered for placement until all three forms are returned to AMBUCS: this form filled out by the recipient/guardian, as well as the Assessment Form and Tryke Selection Form filled out by the therapist.


Amtryke Assessment Form

(Must be filled out completely by therapist)

Recipient's Name: _____
 Age: _____ Weight (lbs.):* _____ Height (inches):* _____
 Diagnosis(es):* _____

*This private information is only used to help appropriately fit the rider.

RIDER'S MEASUREMENTS				
Arm Measurements (inches)				Total Length
Left	A to B:		B to C:	
Right	A to B:		B to C:	
Trunk	A to D:			
Leg Measurements (inches)				Total Length
Left	D to E:		E to F:	
Right	D to E:		E to F:	



A	Center of Shoulder
B	Center of Elbow
C	Center of Digit Crease
D	Center of Hip
E	Center of Knee
F	Bottom of Foot

Arm Length & Leg Length Measurements are critical to correct Amtryke Selection
 Sizing Chart is available online:
www.ambucs.org/riders/wish-list/sizing-chart/

Notes on Provided Measurements (if any):

Helmet Sizing	
Size	Measurement (head circumference)
XXS	18.5" to 19.5"
XS	20.5" to 22"
S/M	22" to 23.6"
L/XL	23.6" to 25.75"

Therapist Name: _____ Is this the treating therapist? Yes No
 Credentials: _____
 Phone: _____ Email: _____
 Facility Name: _____
 Street Address: _____ City: _____ State: _____ Zip: _____
 Is this facility an Amtryke Evaluation Site? Yes No Not Sure

Therapist comments concerning recipient or goals:

This request/assessment is directed to:
 Local AMBUCS Chapter Name: _____
 National Wish List (AMBUCS Resource Center)

By signing below, you are signifying that in your professional opinion this rider would benefit from an Amtryke. You assume no liability.

Therapist Signature: _____ **Date:** _____

Ship Amtryke To			
Name/Facility: _____			Phone: _____
Street Address: _____			
City: _____	State: _____	Zip: _____	

Recipient's Name:

Thanks for choosing an Amtryke therapeutic tricycle!

In order to accommodate the widest variety of people and disabilities, Amtryke offers many tryke models and each can be customized in a variety of ways. Following the steps below will help you choose to the perfect tryke for your client from what might seem like a dizzying array of options.

Remember you can always refer to our website, www.amtrykestore.org, or the Amtryke catalogue for more information and product images.

Step 1: Fill out the first page of the Amtryke Assessment Form.

Step 2: Choose the way the tryke will be propelled: **Hand & Foot, Foot, or Hand**. Your choice should be based on the rider's ability and therapy goals.

Hand & Foot trykes improve coordination, strength and range of motion. Using all four extremities helps with weakness in any area, even general weakness, and can positively affect tone.

Foot trykes were developed in response to requests from therapists for a traditional tricycle for riders with special needs.

Hand trykes are designed for persons whose lower limbs lack function or those who need special therapy for the upper extremity.

Step 3: Take rider's measurements from the front of this Assessment Form and compare them to our Sizing Chart. (This will narrow the choices considerably)

Step 4: Choose a drive. (If it doesn't mention a choice, then ignore this step.)

Amtrykes come with two drive possibilities: **fixed** drive or **geared** drive. Tryke models have been carefully designed so the drive matches the therapeutic goals of the equipment; therefore all drives are not available on all trykes.

A **fixed** drive, commonly known as a 'fixie,' works on a mechanical level to help individuals make a full pedal rotation. The foot crank is constantly in motion for full therapeutic effect. Coasting is not possible with a positive drive; when limb motion stops, the bike does as well.

Riders of **geared** drive trykes must be able to make a complete pedal rotation on their own. This tryke is suitable for riders who need help with balance and a stable sitting position and who have the cognitive ability to steer successfully and apply the coaster brake or hand brake. The key feature of geared trykes is the ability to coast.

Step 5: Chose any adaptations and/or accessories needed by the rider. *Each tryke model can only be customized in the ways noted in its own model section in the Tryke Selection Forms or with the generic accessories listed below. If a customization option is not listed, it is because of design or other practical constraints in offering it on a particular model.*

GENERIC ACCESSORIES (not model specific)

Safety Items	<input type="checkbox"/> Safety Flag	<input type="checkbox"/> Water Bottle w Cage	<input type="checkbox"/> License Plate
Leg & Foot Items	<input type="checkbox"/> Foot Cups (pair): <input type="checkbox"/> Small <input type="checkbox"/> Medium	<input type="checkbox"/> Knee Adductor Strap: <input type="checkbox"/> Small <input type="checkbox"/> Large <input type="checkbox"/> Medium	<input type="checkbox"/> Pedal Block (1 = ¾") _____ qty
Hand Items	<input type="checkbox"/> Variable Range of Motion Kit <i>(only for Hand & Foot Cycles)</i>	<input type="checkbox"/> Wrist Wraps (Includes right & left) <input type="checkbox"/> X-Small <input type="checkbox"/> Large <input type="checkbox"/> Small <input type="checkbox"/> Medium	<input type="checkbox"/> Wrist Brace Mitt: <input type="checkbox"/> Right <input type="checkbox"/> Left <input type="checkbox"/> XX-Small <input type="checkbox"/> Medium <input type="checkbox"/> X-Small <input type="checkbox"/> Small

Recipient's Name: _____

FOOT TRYKES

<input type="checkbox"/> 1410 FOOT TRYKE – (10" wheels, fixed drive) 50-FC-0100 Arm Length 12-20", Leg Length 15-21"			
<input type="checkbox"/> Option #1 – Standard Seating System: Blue Bucket Seat			
<input type="checkbox"/> Option #2 – Alternate Seating System: Snappy Seat System			
Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Laterals <input type="checkbox"/> Head Rest			
<input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back** <input type="checkbox"/> Separator Cube			
<input type="checkbox"/> Option #3 – Alternate Seating System: OLD Blue Bucket Seat** <i>(Less expensive option for very small riders. Cannot be used with an H-Harness.)</i>			
1410 Accessories: <input type="checkbox"/> Rear Steering Kit <input type="checkbox"/> Pull/Steering Bar <input type="checkbox"/> H-Harness <input type="checkbox"/> 8" <input type="checkbox"/> 11.5"			
<input type="checkbox"/> Push Bar <input type="checkbox"/> Heavy-Duty Push Bar			

<input type="checkbox"/> 1412 ProSeries FOOT TRYKE (12" wheels, fixed drive) – 50-FC-1412 Arm Length 14-22", Leg Length 21-26"			
<input type="checkbox"/> Option #1 – Standard Seating System: Saddle Seat/1600 Simple Seatback			
Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Skinny Saddle Seat (8.5x6")			
<input type="checkbox"/> Option #2 – Alternate Seating System: 1400 ProSeries Seat Back System (push grip, backpad, saddle seat & 2 laterals)			
Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Skinny Saddle Seat (8.5x6")			
Trunk Support Options: <input type="checkbox"/> ProSeries Full Padded Back <input type="checkbox"/> ProSeries Head Rest <input type="checkbox"/> ProSeries Lumbar Pad			
<input type="checkbox"/> Recumbent 10 Degree Seat Post			
<input type="checkbox"/> Option #3 – Alternate Seating System: Gray Bucket Seat			
<input type="checkbox"/> Option #4 – Alternate Seating System: Snappy Seat System <i>(too small for many 1412 riders; check seat measurements.)</i>			
Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Laterals <input type="checkbox"/> Head Rest			
<input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back** <input type="checkbox"/> Separator Cube			
1412 Accessories: <input type="checkbox"/> Pull/Steering Bar <input type="checkbox"/> Toe Pedal Leveler Pulley <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5"			
<input type="checkbox"/> Rear Steering Kit <input type="checkbox"/> Front Hand Brake Kit (cannot be used with Rear Steering)			
<input type="checkbox"/> Push Bar <input type="checkbox"/> Hi-Rise Handlebars: <input type="checkbox"/> 5" (BMX) <input type="checkbox"/> 8" <input type="checkbox"/> Exercise Pedals			
<input type="checkbox"/> Heavy-Duty Push Bar <input type="checkbox"/> Knee Separator: <input type="checkbox"/> 2" <input type="checkbox"/> 5" <input type="checkbox"/> 9" Extender Tube			

*cannot be used in conjunction with rear steering. **cannot be used in conjunction with an h-harness, head rest or laterals.

<input type="checkbox"/> NEW! High Functioning Set Up (Both Options Allow Coasting) Can be used with 1416, 1420, 1420XL (Rider must have the ability to steer, apply brake and make a complete pedal rotation on their own.) Comes Standard with Saddle Seat, 1600 Simple Seat Back, Hi-Rise Handlebars and Exercise pedals. MUST choose Option 1 or 2.			
Choose Tryke: <input type="checkbox"/> 1416 <input type="checkbox"/> 1420 <input type="checkbox"/> 1420XL			
<input type="checkbox"/> Option #1: (switch to existing freewheel hub and add Disk Brake Kit - <i>necessary</i> for safety) Adds the ability to coast			
<input type="checkbox"/> Option #2: (3-Speed Kit and coaster brake) Adds the ability to coast and switch gears for hilly terrain – <i>Replaces the 1600 Series</i>			
Standard Seating System (for either option): Saddle Seat/1600 Simple Seat Back			
Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Tractor Seat with bracket <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Med (11.5x10") <input type="checkbox"/> Large(14x13")			
<input type="checkbox"/> Medium Saddle (10.5x10") <input type="checkbox"/> Large Saddle (13x12")			
Adaptive Accessories: <input type="checkbox"/> Expanding Pedals <input type="checkbox"/> Telescoping Loop Handlebars <input type="checkbox"/> 1400 Seat Back Set (w/ laterals)			

Recipient's Name:

FOOT TRYKES

<input type="checkbox"/> 1416 ProSeries FOOT TRYKE (16" wheels, fixed or freewheel drive) – 50-FC-1416 Arm Length 14-22", Leg Length 24-32"
<input type="checkbox"/> Option #1 – Standard Seating System: 1400 ProSeries Seat Back System (push grip, backpad, saddle seat & 2 laterals) Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Tractor Seat with bracket <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Med (11.5x10") <input type="checkbox"/> Large(14x13") Trunk Support Options: <input type="checkbox"/> ProSeries Full Padded Back <input type="checkbox"/> ProSeries Head Rest <input type="checkbox"/> ProSeries Lumbar Pad <input type="checkbox"/> Recumbent 10 Degree Seat Post
<input type="checkbox"/> Option #2 – Alternate Seating System: Saddle Seat/1600 Simple Seat Back Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Tractor Seat with bracket <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Med (11.5x10") <input type="checkbox"/> Large(14x13")
<input type="checkbox"/> Option #3 – Alternate Seating System: Gray Bucket Seat
1416 Accessories: <input type="checkbox"/> Pull/Steering Bar <input type="checkbox"/> Toe Pedal Leveler Pulley <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5" <input type="checkbox"/> Rear Steering Kit <input type="checkbox"/> Calf & Leg Supports <input type="checkbox"/> Exercise Pedals <input type="checkbox"/> Push Bar <input type="checkbox"/> Hi-Rise Handlebars: <input type="checkbox"/> 8" <input type="checkbox"/> 12" <input type="checkbox"/> Heavy-Duty Push Bar <input type="checkbox"/> Knee Separator: <input type="checkbox"/> 2" <input type="checkbox"/> 5" <input type="checkbox"/> 9" extender tube

<input type="checkbox"/> 1420 ProSeries FOOT TRYKE (20" wheels, fixed or freewheel drive) – 50-FC-1420 Arm Length 18-28", Leg Length 27-35"
<input type="checkbox"/> Option #1 – Standard Seating System: 1400 ProSeries Seat Back System (push grip, backpad, saddle seat & 2 laterals) Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Tractor Seat w/ bracket <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Med (11.5x10") <input type="checkbox"/> Large(14x13") Trunk Support Options: <input type="checkbox"/> ProSeries Head Rest <input type="checkbox"/> ProSeries Full Padded Back <input type="checkbox"/> ProSeries Lumbar Pad <input type="checkbox"/> Recumbent 10 Degree Seat Post
<input type="checkbox"/> Option #2 – Alternate Seating System: Saddle Seat/1600 Simple Seatback Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Tractor Seat w/ bracket <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Med 11.5x10" <input type="checkbox"/> Large(14x13")
<input type="checkbox"/> Option #3 – Alternate Seating System: Gray Bucket Seat
1420 Accessories: <input type="checkbox"/> Pull/Steering Bar <input type="checkbox"/> Toe Pedal Pulley <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5" <input type="checkbox"/> Rear Steering Kit <input type="checkbox"/> Calf & Leg Supports <input type="checkbox"/> Exercise Pedals <input type="checkbox"/> Push Bar <input type="checkbox"/> 12" Hi-Rise Handlebars <input type="checkbox"/> Heavy-Duty Push Bar <input type="checkbox"/> Knee Separator: <input type="checkbox"/> 2" <input type="checkbox"/> 5" <input type="checkbox"/> 9" extender tube

<input type="checkbox"/> 1420XL ProSeries FOOT TRYKE (20" wheels, fixed or freewheel drive) – 50-FC-1420-XL Arm Length 20-32", Leg Length 28-43"
<input type="checkbox"/> Option #1 – Standard Seating System: 1400 ProSeries Seat Back System (push grip, backpad, saddle seat & 2 laterals) Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Tractor Seat w/ bracket <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Med (11.5x10") <input type="checkbox"/> Large(14x13") Trunk Support Options: <input type="checkbox"/> ProSeries Head Rest <input type="checkbox"/> ProSeries Full Padded Back <input type="checkbox"/> ProSeries Lumbar Pad <input type="checkbox"/> Recumbent 10 Degree Seat Post
<input type="checkbox"/> Option #2 – Alternate Seating System: Saddle Seat/1600 Simple Seat Back Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Tractor Seat with bracket <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Med (11.5x10") <input type="checkbox"/> Large(14x13")
<input type="checkbox"/> Option #3 – Alternate Seating System: Wheelchair Seat Wheelchair Seat Accessories: <input type="checkbox"/> Wheelchair Seat Bar Ends <input type="checkbox"/> Swing Away Arms
1420 XL Accessories: <input type="checkbox"/> Pull/Steering Bar <input type="checkbox"/> Toe Pedal Pulley <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5" <input type="checkbox"/> Rear Steering Kit <input type="checkbox"/> Calf & Leg Supports <input type="checkbox"/> Exercise Pedals <input type="checkbox"/> Push Bar <input type="checkbox"/> Hi-Rise Handlebars <input type="checkbox"/> 12" <input type="checkbox"/> 23" <input type="checkbox"/> Heavy-Duty Push Bar <input type="checkbox"/> Knee Separator: <input type="checkbox"/> 2" <input type="checkbox"/> 5" <input type="checkbox"/> 9" extender tube

Recipient's Name: _____

FOOT TRYKES

2722 FOOT TRYKE (with 24" wheels, single speed, hand & coaster brake) – **50-FC-2722** Arm Length 22-30", Leg Length 28-36"

Advanced Riders Only; must be able to make a pedal revolution on their own

Option #1 – Standard Seating System: Large Saddle Seat (13x12")

Seat Bottom Alternates: Bench Seat Med Saddle (10.5x10) Tractor Seat w/Bracket: 19" 23"
 Pommel Saddle Seat Medium (11.5x10") Large(14x13")

Option #2 – Alternate Seating System: 1400 ProSeries Seat Back (push grip, backpad, saddle seat (13x12") & 2 laterals)

Seat Bottom Options: Bench Seat Med Saddle (10.5x10) Tractor Seat w/Bracket: 19" 23"
 Pommel Saddle Seat Medium (11.5x10") Large(14x13")

Option #3 – Alternate Seating System: Large Saddle Seat (13x12")/ 1600 Simple Seat Back

Seat Bottom Options: Bench Seat Med Saddle (10.5x10) Tractor Seat w/Bracket: 19" 23"
 Pommel Saddle Seat Medium (11.5x10") Large(14x13")

2722 Accessories: Exercise Pedals Expanding Pedals Heel Pedal Leveler Pulley (only w/ Expanding Pedals)
 23" Hi-Rise Handlebars Calf & Leg Supports (only w/ Expanding Pedals)
 ProSeries Handlebars Stem Riser 7-Speed Kit

JT-2000 FOOT TRYKE (Recumbent 14-Speed) – **50-FC-2000** (advanced riders only) Arm Length 20-28", Leg Length 30-41"

JT-2000 Accessories: Exercise Pedals XL Exercise Pedals Expanding Pedals Digital Speedometer
 Dual Hand Brake Toe Clips Basket Rearview Mirror

JT-2300-USS FOOT TRYKE (Recumbent 14-Speed with **Under Seat Steering**) – **50-FC-2300-USS** (advanced riders only)

Arm Length 20-28", Leg Length 30-41"

JT-2300-USS Accessories: Exercise Pedals XL Exercise Pedals Expanding Pedals Digital Speedometer
 Dual Hand Brake Toe Clips Basket Rearview Mirror

TP-3000 TADPOLE FOOT TRYKE (Recumbent 21-Speed) – **50-FC-3000** (advanced riders only) Arm Length 17-29",

Leg Length 36-45"

Tadpole Accessories: Exercise Pedals XL Exercise Pedals Expanding Pedals Digital Speedometer
 Dual Hand Brake Toe Clips Rearview Mirror

Therapist Assembly Notes or Comments (if any):

Therapist Assembly Notes or Comments (if any):

Recipient's Name: _____

HAND-FOOT TRYKES

<input type="checkbox"/> AM-10 HAND-FOOT TRYKE (10" wheels & 2.5" crank arms) 50-HFC-0105 Arm length 13-17", leg length 15-21"
<input type="checkbox"/> Option #1 – Standard Seating System: Blue Bucket Seat
<input type="checkbox"/> Option #2 – Upgraded Seating System: Snappy Seat System Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Laterals <input type="checkbox"/> Head Rest <input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back**
AM-10 Accessories: <input type="checkbox"/> Rear Steering Kit <input type="checkbox"/> Pull/Steering Bar <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5" <input type="checkbox"/> Push Bar <input type="checkbox"/> Heavy-Duty Push Bar <input type="checkbox"/> Cube Separator

<input type="checkbox"/> AM-12Small HAND-FOOT TRYKE (12" wheels & 3" crank arms) 50-HFC-0110 Arm length 15-20", leg length 19-24" Tryke comes standard with rear steering!
<input type="checkbox"/> Option #1 – Standard Seating System: Bench Seat & 1600 Simple Seat Back
<input type="checkbox"/> Option #2 – Upgraded Seating System: Snappy Seat System Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Laterals <input type="checkbox"/> Head Rest <input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back**
<input type="checkbox"/> Option #3 – Upgraded Seating System: Blue Bucket Seat
AM-12S Accessories: <input type="checkbox"/> Pull/Steering Bar <input type="checkbox"/> Toe Pedal Leveler Pulley <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5" <input type="checkbox"/> Cube Separator <input type="checkbox"/> Exercise pedals

<input type="checkbox"/> AM-12 HAND-FOOT TRYKE (12" wheels & 4" crank arms) 50-HFC-0210 Arm length 14-23", leg length 21-29" Tryke comes standard with rear steering!
<input type="checkbox"/> Option #1 – Standard Seating System: Saddle Seat/1600 Simple Seat Back Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Large(14x13")
<input type="checkbox"/> Option #2 – Upgraded Seating System: ProSeries Seat Back System (push grip, backpad, saddle seat & 2 laterals) Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Large(14x13") Trunk Support Upgrades: <input type="checkbox"/> ProSeries Full Padded Back <input type="checkbox"/> ProSeries Head Rest <input type="checkbox"/> ProSeries Lumbar Pad
<input type="checkbox"/> Option #3 – Upgraded Seating System: Gray Bucket Seat
<input type="checkbox"/> Option #4 – Upgraded Seating System: Snappy Seat System Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Laterals <input type="checkbox"/> Head Rest <input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back** <input type="checkbox"/> Cube Separator
AM-12 Accessories: <input type="checkbox"/> Pull/Steering Bar <input type="checkbox"/> Toe Pedal Leveler Pulley <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5" <input type="checkbox"/> Knee Separator: <input type="checkbox"/> 2" <input type="checkbox"/> 5" <input type="checkbox"/> Exercise pedals

<input type="checkbox"/> AM-16 HAND-FOOT TRYKE (16" wheels & 5" crank arms) 50-HFC-0411 Arm length 18-27", leg length 24-36" Now comes standard with rear steering!
<input type="checkbox"/> Option #1 – Standard Seating System: Saddle Seat/1600 Simple Seat Back Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Tractor Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Large(14x13")
<input type="checkbox"/> Option #2 – Upgraded Seating System: 1400 ProSeries Seat Back System (push grip, backpad, saddle seat & 2 laterals) Seat Bottom Upgrades: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Tractor Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Large(14x13") Trunk Support Upgrades: <input type="checkbox"/> ProSeries Full Padded Back <input type="checkbox"/> ProSeries Head Rest <input type="checkbox"/> ProSeries Lumbar Pad
<input type="checkbox"/> Option #3 – Upgraded Seating System: Gray Bucket Seat
AM-16 Accessories: <input type="checkbox"/> Pull/Steering Bar <input type="checkbox"/> Toe Pedal Leveler Pulley <input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5" <input type="checkbox"/> Knee Separator: <input type="checkbox"/> 2" <input type="checkbox"/> 5" <input type="checkbox"/> 9" Extender Tube <input type="checkbox"/> Exercise pedals <input type="checkbox"/> XL Exercise pedals

*cannot be used in conjunction with rear steering. **cannot be used in conjunction with an h-harness, head rest or laterals.

Please mail, email or fax completed form to your local chapter or AMBUCS Resource Center

Mail: P.O. Box 5127, High Point, NC 27262 Email: wishlist@ambucs.org Fax: 336.852.6830

This form, Request & Liability Waiver Form and Assessment Form must be received before placement is considered.

(800) 838-1845

ambucs.org

Recipient's Name:

New! AM-20 HAND-FOOT TRYKE (Recumbent Single Speed) – **50-HFC-0610** Arm length 18-30", leg length 32-46"

Tryke is currently being used for research!

AM-20 Accessories: Exercise Pedals XL Exercise Pedals Expanding Pedals
 Wheelchair Seat Bar Ends Swing Away Arms Calf & Leg Supports (Expanding pedals)

Therapist Assembly Notes or Comments (if any):

Recipient's Name:	
--------------------------	--

HAND TRYKES

<input type="checkbox"/> AM-10 HAND TRYKE (12" wheels, 2.5" crank arms & permanent foot platform) 50-HFC-0105 Arm length 13-17"			
<input type="checkbox"/> Option #1 – Standard Seating System: Blue Bucket Seat			
<input type="checkbox"/> Option #2 – Upgraded Seating System: Snappy Seat System			
Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Laterals <input type="checkbox"/> Head Rest			
<input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back**			
AM-10 Accessories:	<input type="checkbox"/> Rear Steering Kit	<input type="checkbox"/> Pull/Steering Bar	<input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5"
	<input type="checkbox"/> Push Bar	<input type="checkbox"/> Heavy-Duty Push Bar	<input type="checkbox"/> Separator Cube

<input type="checkbox"/> AM-12Small HAND TRYKE (12" wheels, 3" crank arms & permanent foot platform) 50-HFC-0110 Arm length 15-20"			
Tryke comes standard with rear steering!			
<input type="checkbox"/> Option #1 – Standard Seating System: Bench Seat & 1600 Simple Seat Back			
<input type="checkbox"/> Option #2 – Upgraded Seating System: Snappy Seat System			
Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Laterals <input type="checkbox"/> Head Rest			
<input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back**			
<input type="checkbox"/> Option #3 – Alternate Seating System: Blue Bucket Seat			
AM-12S Accessories:	<input type="checkbox"/> Pull/Steering Bar	<input type="checkbox"/> Toe Pedal Leveler Pulley	<input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5"
			<input type="checkbox"/> Separator Cube

<input type="checkbox"/> AM-12 HAND TRYKE – (12" wheels, 4" crank arms & permanent foot platform) 50-HC-0210 Arm Length 14-23"			
Tryke comes standard with rear steering!			
<input type="checkbox"/> Option #1 – Standard Seating System: Saddle Seat/1600 Simple Seat Back			
Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Large(14x13")			
<input type="checkbox"/> Option #2 – Alternate Seating System: 1400 ProSeries Seat Back System (push grip, backpad, saddle seat & 2 laterals)			
Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Large(14x13")			
Trunk Support Options: <input type="checkbox"/> ProSeries Full Padded Back <input type="checkbox"/> ProSeries Head Rest <input type="checkbox"/> ProSeries Lumbar Pad			
<input type="checkbox"/> Option #3 – Alternate Seating System: Gray Bucket Seat			
<input type="checkbox"/> Option #4 – Alternate Seating System: Snappy Seat System			
Snappy Seat Accessories: <input type="checkbox"/> Notched Seat Bottom Cushion* <input type="checkbox"/> Laterals <input type="checkbox"/> Head Rest			
<input type="checkbox"/> Notched Wide Bottom Cushion* <input type="checkbox"/> XL Seat Back** <input type="checkbox"/> Separator Cube			
AM-12 Accessories:	<input type="checkbox"/> Pull/Steering Bar	<input type="checkbox"/> Toe Pedal Leveler Pulley	<input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5"
	<input type="checkbox"/> Knee Separator: <input type="checkbox"/> 2" <input type="checkbox"/> 5"	<input type="checkbox"/> 9" Extender Tube	

*cannot be used in conjunction with rear steering. **cannot be used in conjunction with an h-harness, head rest or laterals.

<input type="checkbox"/> AM-16 HAND TRYKE (16" wheels, 5" crank arms & permanent foot platform) 50-HFC-0411 Arm length 18-27"			
Tryke comes standard with rear steering!			
<input type="checkbox"/> Option #1 – Standard Seating System: Saddle Seat/1600 Simple Seat Back			
Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Large(14x13")			
<input type="checkbox"/> Option #2 – Upgraded Seating System: 1400 ProSeries Seat Back System (push grip, backpad, saddle seat & 2 laterals)			
Seat Bottom Alternates: <input type="checkbox"/> Bench Seat <input type="checkbox"/> Padded Tractor Seat <input type="checkbox"/> Pommel Saddle Seat <input type="checkbox"/> Medium (11.5x10") <input type="checkbox"/> Large(14x13")			
Trunk Support Options: <input type="checkbox"/> ProSeries Full Padded Back <input type="checkbox"/> ProSeries Head Rest <input type="checkbox"/> ProSeries Lumbar Pad			
<input type="checkbox"/> Option #3 – Upgraded Seating System: Gray Bucket Seat			
AM-16 Accessories:	<input type="checkbox"/> Pull/Steering Bar	<input type="checkbox"/> Toe Pedal Leveler Pulley	<input type="checkbox"/> H-Harness: <input type="checkbox"/> 8" <input type="checkbox"/> 11.5"
	<input type="checkbox"/> Knee Separator: <input type="checkbox"/> 2" <input type="checkbox"/> 5"	<input type="checkbox"/> 9" Extender Tube	

Please mail, email or fax completed form to your local chapter or the AMBUCS Resource Center
 Mail: P.O. Box 5127, High Point, NC 27262 Email: wishlist@ambucs.org Fax: 336.852.6830

(800) 838-1845
ambucs.org

This form, Request Form & Liability Waiver Form and Assessment Form must be received before Amtryke placement is considered.

Recipient's Name: _____

HAND TRYKES

1024 HAND TRYKE *fka 1424* – (24" wheels & 3-Speed drive train) **50-HC-1024** Arm length 22-26"

1024 Accessories: Wheelchair Seat Bar Ends Swing Away Arms Remote location hand brake & shifting kit

HP-1000 RECREATION HAND TRYKE – (16" front wheel & 7-Speed grip shift) **50-HC-1000** Arm length 22-26"


HP-1000 Accessories: Wheelchair Seat Bar Ends Swing Away Arms

Therapist Assembly Notes or Comments (if any):

Amtryke Sizing Chart

TRYKE TYPE (How will the tryke be propelled?)	RIDER LEG LENGTH (Inches from center of hip to bottom of shoe.)	RIDER ARM LENGTH (Inches from middle of shoulder to center of digit crease.)	MODEL	RIDER WEIGHT (Pounds)	RIDER MAX HEIGHT (Inches)	TRYKE WEIGHT (Pounds)	WHEEL SIZE (Inches)	TTRYKE HEIGHT (inches)	TRYKE LEGNTH (inches)	TRYKE WIDTH (Inches)
Hand & Foot	15-21	13-17	AM-10	55	40	45	10	24	38	21
	19-24	15-20	AM-12S	150	40	45	12	27	38	24
	21-29	14-23	AM-12	150	47	45	12	36	60	32
	24-36	18-27	AM-16	175	66	55	16	36	68	33
	32-46	18-30	AM-20	250	74	80	20	45	72	32
Foot	15-21	12-20	1410	55	40	45	10	24	38	21
	21-26	14-22	1412	125	42	72	12	13	43	27
	24-30	14-22	1416	175	60	74	16	49	58	30
	27-33	18-28	1420	250	68	74	20	50	64	30
	28-43	20-32	1420XL	275	76	89	20	43	72	29
	28-36	22-30	2722	275	72	76	24	42	60	30
	30-41	20-28	JT-2000/ JT2300USS	250	74	80	20	48	72	32
	36-45	17-29	TP-3000	300	75	47	20	33	63-70	33
Hand	22-41	22-26	1024	250	72	74	20	45	75	32
	22-41	22-26	HP-1000	250	74	80	Front 16 Rear 20	45	72	32
	All trykes in the Hand & Foot section can be converted to Hand trykes – except the AM-20.									

- A Center of Shoulder
- B Center of Elbow
- C Center of Digit Crease
- D Center of Hip (greater trochanter)
- E Center of Knee
- F Bottom of Foot



RIDER'S MEASUREMENTS	
Arm Measurements (inches) Total Length	
Left	A to B + B to C = _____
Right	A to B + B to C = _____
Trunk	A to D = _____
Leg Measurements (inches) Total Length	
Left	D to E + E to F = _____
Right	D to E + E to F = _____

HELMET SIZING	
Youth Sizes	Head Circumference Inches
XXS	18.5 to 19.5
XS	20.5 to 22
S/M	22 to 23.6
L/XL	23.6 to 25.75
Adult Sizes	Head Circumference Inches
S/M	22 to 23.6
L/XL	23.6 to 25.75